

LYMPH 0 FIT™



Lymphatic relief arm sleeve Halifax 1115

Indication


- Lymphoedema in the arm
- After arm lift surgery
- After liposuction

Application

- The special 3D dot texture creates a massaging and draining effect on the dermis and the smaller lymph vessels
- Suitable for prophylactic use with all lymph issues

The arm sleeve is easy to put on

Flat seams with no pressure points – helps to prevent lymphoedema



Special cuff on the sleeve end to prevent slipping and pressure points

Open edging (clean cut) on the fingers, with no thumb section, to avoid lymphatic congestion

Tip: If you are between two sizes, always opt for the larger size.

SIZE	CIRCUMFERENCE WRIST
0	12 – 13 cm
1	14 – 15 cm
2	16 – 17 cm
3	17 – 18 cm
4	19 – 20 cm
5	20 – 21 cm

SIZE	CIRCUMFERENCE BICEPS
0	22 – 27 cm
1	24 – 29 cm
2	26 – 32 cm
3	28 – 34 cm
4	30 – 36 cm
5	32 – 38 cm

SIZE	CIRCUMFERENCE ELBOW
0	22 – 23 cm
1	23 – 25 cm
2	24 – 27 cm
3	25 – 29 cm
4	26 – 31 cm
5	27 – 33 cm

AVAILABLE COLOURS

001	black
007	sand



FABRIC

71% PA, 29% EL

40° C washable, wash on delicates cycle